

## ChIP-Seq Library Prep

## Version 1.0

## Date:

## Starting Materials

- 1.
- 2.
- 3.
- 4.
- 5.
- 6.
- 7.
- 8.
- 9.
- 10.

## End Repair (NEB Next Kit)

30ul DNA

5ul Neb Next Buffer

2.5ul Enzyme Mix

12.5ul H2O

50ul Total

Incubate 30 minutes at 20 degrees in PCR machine

Purify with 1.8 volumes ampure beads (90ul) - elute in 43ul H2O

## A-tail

42ul DNA

5ul NEB #2

1ul dATP

2ul Klenow exo-

50uL Total

Incubate 30 minutes at 37 degrees

Purify with 1.8 volumes ampure beads (90uL) - elute 35ul H2O

### Adaptor Ligation (NEBNext Kit)

34ul DNA

1ul adapter mix (1uM if <10ng starting material or 1.5uM if <50ng)

10ul 5X ligase buffer

5ul Ligase

50ul Total

Incubate 30 minutes at 20 degrees

Purify 2X with ampure beads

*1st round - 1 volume ampure; elute in 50ul*

2nd round - 0.8 volume ampure; elute in 51ul

**Check Linearity with qPCR - use 1/20 of sample with MC703-MC704**

Use threshold value at midpt of amplification to calculate final number of cycles to use

—

PCR Enrichment with NEB Next High Fidelity

DNA

2X NebNext High Fidelity Mix 25ul

20uM PE primer 1(MC703) 1ul

20uM PE primer 2(MC704) 1ul

H2O

**Total 50ul**

Cycling Parameters

98X 30"

98X 10"

65X 30"

72X 30"

72X 5'

Hold at 4 degrees

Purify Samples with 50ul Ampure Beads and elute with 20ul H2O

Run 1ul of this on 2% agarose gel - stained with sybr gold (30 minutes)

Spec 1ul on Qubit